# Genebank - Cryopreservation of Potato - OP106

**Date** : September, 2015

**Reference** : Vollmer, R.; Villagaray, R.; Egusquiza, V.; Panta, A.; Barkley, N.; Ellis, D. 2016.  
Genebank: Cryopreservation of potato, v.2. CIP-OP106. 30 pps. Not for general distribution. For most current version of the procedure, please contact Rainer Vollmer (r.vollmer@cgiar.org) 1-34 pp.

Not for general distribution.

`$metadata.from("CONTACT")`

Genebank - Cryopreservation of Potato, v.2. CIP-OP106 by International Potato Center (CIP) is licensed under a Creative Commons Attribution-ShareAlike 4.0 International License. Permissions beyond the scope of this license may be available at <http://www.cipotato.org/>

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<b>TITLE</b>	<b>Genebank - Cryopreservation of Potato - OP106</b>
<b>OWNER *</b>	Head Cryopreservation Unit
<b>APPROVER *</b>	Head Genebank
<b>APPROVAL DATE *</b>	<b>Pending approval, under revision and edition.</b>
<b>LAST REVIEW DATE *</b>	September, 2015
<b>ISSUE DATE</b>	Apr 26, 2016 14:52
<b>CONTRIBUTORS *</b>	Vollmer, Rainer (CIP), Villagaray, Rosalva (CIP), Egusquiza, Veronica (CIP), Panta, Ana (CIP), Ellis, David (CIP)
<b>CITATION</b>	106
<b>KEYWORDS</b>	procedure, non-accredited
<b>DOCUMENT ID *</b>	OP106
<b>VERSION NUMBER</b>	V2
<b>COMPETENT PERSONNEL *</b>	Villagaray, Rosalva (CIP), Egusquiza, Veronica (CIP), Cardenas, Jose (CIP), Castro, Mario (CIP), Torres, Maria, Panta, Ana (CIP), Vollmer, Rainer (CIP)
<b>COMPETENCE TEST *</b>	
<b>DATE OF NEXT COMPETENCE TEST *</b>	
<b>CITATION OF THIS DOCUMENT *</b>	Vollmer, R.; Villagaray, R.; Egusquiza, V.; Panta, A.; Barkley, N.; Ellis, D. 2016. Genebank: Cryopreservation of potato, v.2. CIP-OP106. 30 pps. Not for general distribution. For most current version of the procedure, please contact Rainer Vollmer (r.vollmer@cgiar.org)

- All fields marked with \* are required

## ENGLISH VERSION

## INTRODUCTION

---

Cryopreservation is an effective strategy for the long-term conservation of plant genetic resources collections. Plant material is stored at ultra low temperatures in liquid nitrogen ( $-196^{\circ}\text{C}$ ) and at such low temperatures physiological and chemical activities are reduced to a point where germplasm should remain viable for centuries. As of April 2016, CIP's cryobank holds about 1450 potato accessions. By the end of 2016, the potato cryobank will be a safety back-up for 35% of CIP's native potato varieties.

## SCOPE

---

The following is CIP's standard operating procedure for the cryopreservation of potato using the PVS2-droplet vitrification method.

# HEALTH & SAFETY

1. CIP's Health & Safety Manual
2. CIP's Manual for "Safe Use, Handling and Storage of Liquid Nitrogen"
3. Protocol for the "Supply, Monitoring and Transfer of Liquid Nitrogen (LN)" at CIP's cryobank
4. Wear safety glasses (or goggles) and mask when working in the laminar flow chamber. Wear a long lab coat during all laboratory work.
5. Do not sterilize tools longer than 30 seconds in the glass bead sterilizer to avoid burns on fingers.
6. Point the scalpel edge away from your body to avoid cutting accidents during installation and deinstallation. Installation and desinstallation should be done smoothly without using much force or pressure.
7. Wear hermetic safety goggles and mask during freezing, thawing, and transfer of cryovials. People working near the person who is freezing, thawing, or transferring cryovials also have to wear hermetic safety goggles.
8. Wear cryogenic gloves (lab model), hermetic safety goggles, and cryo-apron when pouring liquid nitrogen (LN) from dewars to freezing containers or pans. People working near the person who is pouring LN also have to wear hermetic safety goggles.
9. Wear cryogenic gloves (industrial), hermetic safety goggles, cryo-apron, and face shield when handling high-pressure tanks or transferring bigger volumes of LN (i.e. filling dewars, refilling cryotanks, dispensing LN with hoses, etc.)
10. Always wear closed-toe shoes, and long pants or skirts which fully cover the legs. Never wear sandals or open-toed shoes in the laboratory or cryobank.
11. Take periodic active breaks when working in a static or fixed position (i.e. do a 2 minute stretching exercise every working hour).
12. When working with a stereoscope, close eyes and focus on different distances every 15 minutes (to reduce eye strain).
13. Wear surgical gloves, mask, and safety glasses when handling chemical products.

## MATERIALS

Equipment	
Barcode printers (desktop and hand-held)	Laminar Flow Chamber
Cryotank for Temporary Storage (Taylor Wharton LS-4800)	Large Cryotank (CHART MVE-819P-190AF-GB equipped with TEC3000) (= <b>cryobank</b> )

<b>Equipment</b>	
Dewar for liquid nitrogen (4 liters)	Medium-sized cryotanks (= working tanks) (CHART MVE XC47/11 and CHART Cryosystem 2000)
Flake Ice Machine	Pocket PC
Laboratory Mixer (Maxi Mix 2)	Stereoscope
<b>Other materials</b>	
Aluminium foil	Glass petri dishes (100 x 15 mm)
Barcode labels for desktop printer (polypropylene, 1.3" x 1.38")	Glass test tubes with caps (25x150 mm)
Barcode labels for hand-held printer (polypropylene, PolyPro 4000D, 2" x 1")	Glass test tubes with caps (sterile) (13x100 mm) (*)
Blades (#10 and #11)	Glass vessel with wide mouth (approx. 250 mL)
Coolbox system (CF 30 Cool Box System)	Long forceps (23 cm, fine point of 1 mm)
Cotton	Magic Touch lab ice pan (capacity: 9 liter)
Cotton forceps (stainless steel)	Parafilm
Cryo Apron	Pencil
Cryobox (for 25 and 100 cryovials)	Pieces of sterile filter papers (two dimensions: 10x10 mm and 25x30 mm) (*)
Cryogenic gloves (Tempshield)	Rubber bulbs (2 mL) for glass pasteur pipette
Cryogenic labels (ZEBRA 8000T)	Safety glasses
Cryogenic marker	Safety goggles (hermetic)
Cryovials (Thermo Scientific, 1.8 mL, round, internal thread, starfoot base) (*)	Sharp-point pliers
Deep petri dishes (100 x 25 mm; disposable) (*)	Small stickers (approx. 1x1 cm)
Digital Chronometer	Standard petri dishes (100x15 mm; disposable)
Dissecting knife handle	Sterile screw test tubes (15 ml, 70 mm, Ø20 mm) (*)
Face shield	Sterile strips of aluminium paper (5x15 mm) (*)
GA7 magenta vessels (sterile; with lid) (*)	Sterile bond paper (*)
Glass pasteur pipettes (5.75 ") (*)	Test tube racks (for screw test tubes, 13x100 mm and 25x150 mm test tubes)

(\*) for details on preparation, usage, sterilization or washing please see the related document RD200.

Chemicals
Flaked ice
Liquid Nitrogen (LN)
Loading Solution (LS-P) (**)
Potato Propagation Medium (MSA) (**)
Plant Vitrification Solution 2 (PVS2-P) (**)
Potato Recovery Medium I (PRM-0.3 M) (**)
Potato Recovery Medium II (PRM-0.2 M) (**)
Potato Recovery Medium III (PRM-0.1 M) (**)
Potato Recovery Medium IV (PRM-25 g) (**)
Rewarming Solution (RS-P) (**)

(\*\*) For composition and preparation of solutions and culture media please see: Genebank - Cryopreservation - Solutions and Culture Media - OP108

## PROCEDURE

### 1. General selection and quality criteria for cryopreservation

**1.1.** Select accessions to be cryopreserved according to the following priorities (in this order): (a) Accessions included in the International Treaty on Plant Genetic Resources for Food and Agriculture (PGRFA), (b) True-to-type identity confirmed, (c) Phytosanitarily clean (HS2) (if no accessions in HS2 are available, select HS2R), (d) Accessions forming part of the core collection.

**1.2.** For each accession and repetition cryopreserve a group of 150 shoot tips. Each vial contains 10 shoot tips. For viability assessment thaw 30 out of 150 shoot tips after a minimum LN exposure of 24 hours. Assess survival and regrowth rates at 30 and 60 days after thawing. Transfer accessions without fungal or bacterial contamination and that meet the minimum viability criteria (points 1.4 to 1.6) to the cryobank and safety back-up tanks. Accept samples as recovered if they grow in an *in vitro* plant with normal attributes, including elongated stem, apex, leaves and root(s).

**1.3.** If the average viability of the three routinely tested vials of the first run is **equal or higher than 30%**, only a single run is required. Place 2 vials in a box for cryo safety back-up and 10 vials into CIP's cryobank. (\*)

**1.4.** If the average viability of the three routinely tested vials of the first run is **between 20% and 30%**, do a second run of 150 shoot tips. The following process will be used based on the results from the second run:

**(a)** If the second run has an average viability of equal or higher than 20%, place both runs in the cryobank and safety back-up tanks. For the first run place 4 vials in a box for cryo safety back-up and 8 vials into CIP's cryobank (\*\*). The distribution of vials between cryobank and safety back-up of the second run depends on its viability rate:

- If the second run has a viability rate of equal or higher than 20% but less than 30%, place 4 vials in a box for cryo safety back-up and 8 vials into CIP's cryobank (\*\*).
- However, if the second run shows a viability of equal or higher than 30%, place 2 vials in a box for cryo safety back-up and 10 vials into CIP's cryobank (\*).

**(b)** If the second run shows an average viability of less than 20%, discard and redo the second run. If the redone second run shows a viability rate of equal or higher than 20%, proceed as described in step 1.4.a. If the redone second run shows a viability of less than 20%, discard the second run, only store the first run, and label the accession with "low viability" in the database.

**1.5.** If the average viability of the routinely tested samples of the first run is **less than 20%**, do not cryopreserve a second run and discard the first run. Label the accession with "low viability" (less than 20% but not 0%) or "cryo recalcitrant" (0% viability) in the database. These accessions will be cryopreserved with an improved protocol in a later time.

If less than 12 vials are stored, since a contamination recheck is required (see step 10), the distribution of vials between cryobank and safety back-up boxes is as follows:

**(\*)**

- 11 vials --> safety back-up box: 2 vials, cryobank box: 9 vials
- 10 vials --> safety back-up box: 2 vials, cryobank box: 8 vials
- 9 vials --> safety back-up box: 2 vials, cryobank box: 7 vials

**(\*\*)**

- 11 vials --> safety back-up box: 3 vials, cryobank box: 8 vials
- 10 vials --> safety back-up box: 3 vials, cryobank box: 7 vials
- 9 vials --> safety back-up box: 3 vials, cryobank box: 6 vials

## **2. In vitro culture of mother plants**

**2.1 First multiplication cycle:** Propagate axillary sections (with 2-3 nodes) of *in vitro* potato plants on MSA medium (see OP 108) placing four explants per 25x150 mm test tube. Propagate a total number of 6-7 test tubes. Propagate with long forceps and sharp No. 10 scalpel on a stack of 2-3 sterile bond papers. Periodically replace the stack to avoid cross contamination. Label each tube with preprinted barcode labels containing the labcode, CIP number, common name of accession, propagation date, health status, batch of culture media, and ID of person who printed the labels. Incubate for 2-3 weeks at 20±2 °C and 85±20 µmol.m-2.s-1 light intensity, provided by cool daylight fluorescent tubes with a 16h/8h photoperiod.

**2.2. Second multiplication cycle:** Propagate axillary sections (with 2-3 nodes) of *in vitro* potato plants on MSA medium (see OP 108) placing four explants per 25x150 mm test tube. Propagate a total number of 14-16 test tubes. Label each tube with preprinted barcode labels containing labcode, CIP number, common name of accession, propagation date, health status, batch of culture media and ID of person

who printed out the labels. Incubate for 3 weeks at  $20\pm 2$  °C and  $85\pm 20$   $\mu\text{mol}\cdot\text{m}^{-2}\cdot\text{s}^{-1}$  light intensity, provided by cool daylight fluorescent tubes with a 16h/8h photoperiod.

**2.3. Third multiplication cycle:** Propagate uni-nodal axillary sections of *in vitro* potato plants on MSA medium (see OP 108), with 25 explants per magenta vessel. Propagate a total number of 8 magenta vessels. Label each magenta vessel with preprinted barcode labels containing labcode, CIP number, common name of accession, propagation date, health status, batch of culture media and ID of person who printed out the labels. Magenta vessels (8) of the same accession are placed on a single GA7 rack to avoid mix-ups and facilitate handling. Incubate for 2 weeks at  $20\pm 2$  °C and  $85\pm 20$   $\mu\text{mol}\cdot\text{m}^{-2}\cdot\text{s}^{-1}$  light intensity, provided by cool daylight fluorescent tubes with a 16h/8h photoperiod.

**2.4 .Final multiplication cycle (= Installation):** Propagate uni-nodal **apical** sections on MSA medium (see OP 108), with 60-70 apical explants per deep petri dish. Propagate a total number of three petri dishes (180-210 explants). Label each petri dish with preprinted barcode labels containing labcode, CIP number, common name of accession, propagation date, health status, batch of culture media and ID of person who printed out the labels. Incubate for 2-3 weeks (genotype dependent) at  $7\pm 1$  °C and  $15\pm 5$   $\mu\text{mol}\cdot\text{m}^{-2}\cdot\text{s}^{-1}$  light intensity, provided by cool daylight fluorescent tubes with a 16h/8h photoperiod.

### **3. Excision of apical shoot tips**

**3.1.** For each accession, distribute 15 small filter papers (10x10 mm) onto a petri dish containing sterile MSA medium (see OP 108).

**3.2.** Under a stereoscope (10x) carefully excise shoot tips from 2 to 3-week-old *in vitro* potato plants. Excise with long forceps and sharp No. 11 on a stack of 2-3 sterile bond papers. Periodically replace the stack to avoid cross contamination. Shoot tips should contain 2-4 leaf primordia, have a length of 0.8-1.2 mm and a width of 0.3-0.7 mm . The reference value of the width of the shoot tip was taken at the height of the meristem, and varies depending on the genotype. For each accession excise 150 shoot tips and place them successively onto the small filter papers until each of the 15 filter papers contains 10 shoot tips (first place one shoot tip on each of the 15 papers, then a second on each, than a third, and so on). At the end of the cryopreservation process, 120 shoot tips will be conserved in liquid nitrogen for long-term conservation and 30 will be thawed to determine viability (after a minimum of 24 hours in liquid nitrogen).

### **4. Preparative steps for cryoprotection**

**4.1.** Disinfect rubber bulbs (2 mL) in a glass vessel with alcohol (96%). Set up sterile standardized glass pasteur pipettes for Loading Solution (LS) and Plant Vitrification Solution 2 (PVS2) with its rubber bulbs and store pipettes individually in sterile 15 ml screw test tubes. Always pipette LS and PVS2 with its corresponding standardized glass pasteur pipette to avoid mixing solutions and problems during shoot tip handling. For composition and preparation of LS and PVS2 please see OP 108. If the pipette tip comes in contact with non-sterile items or absorbs an improper solution, immediately replace the pipette with a new sterile standardized pipette of the same type (LS or PVS2).

**4.2.** Fill lab ice pan (9 liter capacity) with flaked ice. Homogenize LS and PVS2 with a laboratory mixer (Maxi Mix II). Place Eppendorf containing PVS2 on flaked ice (0 °C) at a depth of approx. 3/4 of the

Eppendorf. Print 15 small labels for the LS/PVS2 test tubes (15 mL) and stick them on the caps of 15 sterile test tubes. The LS-PVS2 labels list the labcode of accession, "LS: " and "PVS2: ". Using the corresponding LS pipette dispense 2 mL of LS per screw test tube.

**4.3.** Using the pocket PC, record the following information in the database: ID of run (1st or 2nd), date of excision and freezing, name of person who has excised and frozen the shoot tips, number of cryopreserved shoot tips, abbreviation of protocol, and notes (if necessary). Using the pocket PC print, 15 cryogenic barcode labels per accession.

The barcode labels should contain the following information: CIP number of the accession, labcode of the accession, freezing date, and abbreviation of protocol. Using a cryogenic marker, label 13 sterile cryotubes with their lab-code and freezing date. Label cryotubes with corresponding barcode labels. Locate open cryovials in the coolbox system holder (CF 30 Cool Box System) and label the corresponding rows of the holder with cryogenic labels.

12 of 15 labels identify the vials assigned to the cryobank and safety copy. 1 label identifies the vial that will contain the 3 aluminum strips (30 samples) for routine viability assessment. 1 label identifies the location (rows) of the accession within the freezing container. 1 label is printed as a back-up.

**4.4.** Place previously sterilized aluminum strips (0.5 x 2.0 cm) in a sterile plastic petri dish. Using cotton forceps gently fold the strips on one end (approx. 3 mm from the end of the strip).

## **5. Cryoprotection: Loading (LS) and Plant Vitrification Solution 2 (PVS2)**

**5.1.** Using cotton forceps, introduce the small filter papers holding 10 apical shoot tips into test tubes with LS. Shake smoothly. Verify that all shoot tips are submerged in LS. Treat shoot tips for 20 minutes with LS. Using a pencil, note the time of finishing the LS treatment on the sticker and monitor the exposure time with a digital clock. Make a break of 1-2 minutes between each test tube.

**5.2.** Place a sterile glass vessel in the flow chamber. After 20 minutes of LS exposure, carefully remove the solution with a pasteur pipette such that the shoot tips remain in the test tube. Discard the used LS in the sterile glass vessel. Using a separate pipette, place approx. 2.0 mL of ice-cooled PVS2 (0 °C) into the tube containing the shoot tips. Close tube with a sterile cap. Verify that all shoot tips are submerged in PVS2. Treat shoot tips for 50 minutes with PVS2 (at 0 °C). Using a pencil, note the time of finishing the PVS2 treatment on the sticker. Place the test tube on flaked ice (0 °C) at a depth of approx. 3/4 of the tube. Monitor the PVS2 exposure time with a digital clock.

## **6. Freezing in liquid nitrogen**

**6.1.** Two minutes before the PVS2 treatment is complete, carefully pour LN from a 4 liter dewar into the coolbox system container (wear cryogenic gloves and safety glasses). If necessary, relocate cryovials (without cap) in the supporting holes of the coolbox system holder. The nitrogen level always should be approx. 1.0 cm above the opening of the cryovials. Also place and submerge the cryovial caps in LN for additional sterilization.

**6.2.** Using the PVS2 pasteur pipette, suck up the 10 shoot tips contained in PVS2, lean the pipette's tip on the bottom of the tube and slowly expel the PVS2 from the pipette. The shoot tips should be contained within a very small volume of PVS2 (approx. 10 µL) in the tip of the pipette.

**6.3.** Place the shoot tips, along with the PVS2 droplet on a thin sterile aluminum strip (5 x 20 mm). The aluminum strip should be folded at one end to facilitate manipulation by cotton forceps.

**6.4.** Quickly plunge the aluminum strip into the coolbox system container with liquid nitrogen (fast freezing) and transfer it to the corresponding labelled cryovial. Move the strip only under the surface of liquid nitrogen. Place only one strip (= 10 shoot tips) per vial. Periodically refill the coolbox system container with liquid nitrogen to assure that the cryovials are always submerged. Place the three aluminum strips (= 30 shoot tips) used for viability assessment together in one single cryovial.

**6.5.** When the freezing process of all accessions to be cryopreserved is finished (during the working day), assign 12 cryovials (= 120 shoot tips ) for cryobanking and safety back-up. Using sharp-pointed pliers, tightly close cryovial caps. The three aluminum strips (= 30 shoot tips) for routine viability assessment will be thawed in RS and cultured on recovery media, after a minimum exposure of 24 hours to LN.

**6.6.** Place the cryo-cane based working tank (CHART MVE XC47/11) next to the laminar flow chamber and open the tank. Using an ultra-low temperature resistant pencil, label cryo-canes with the labcode of the accession they will hold. Canes have a capacity of 4 cryovials each. Using an ultra-low temperature resistant pencil, label 3 cryo-canes per accession with their corresponding labcodes and pre-freeze them in LN (in tank holders). Using sharp-pointed pliers, quickly transfer cryovials (12 cryovials per accession) from coolbox system container to the corresponding pre-frozen aluminum cryocane and immediately submerge the cryocane in LN. Expose the cryovials to room temperature for a minimum of time and quickly re-submerge them in LN. The cryo-cane based working tank contains 6 cylindrical holders with a capacity of approx. 45 cryo-canes per holder. The total capacity of the cryo-cane based working tank is 75 accessions.

Locate the cryovial that contains the three aluminum strips (=30 shoot tips) used for routine viability assessment in a separated holder of the working tank. Regularly check the LN level of the cryo-cane based working tank and refill as needed.

## **7. Thawing in Rewarming Solution (RS)**

**7.1.** Disinfect rubber bulbs (2 mL) in a glass vessel with alcohol (96%). Set up sterile standardized glass pasteur pipette for Rewarming Solution (RS) with its rubber bulbs and store the pipette in a sterile 15 mL test tube. For composition and preparation of RS please see OP 108. If the pipette tip comes in contact with non-sterile objects, immediately replace the pipette with a new sterile standardized pipette.

**7.2.** Thaw in groups of 6 accessions. Fill 18 sterile 15 mL test tubes with RS (4 mL of RS per tube) and arrange them in the test tube rack in 6 rows of 3 tubes each. Each row will correspond to an accession. Using a pocket PC, print 4 small labels per accession. The RS labels contain labcode of accession, and "RS:". 3 of 4 labels are used to identify each of the test tubes of the corresponding accession, and the fourth label identifies the petri dish with sterile Potato Recovery Medium I (PRM-0.3 M, see OP108) to which the shoot tips will be transferred to.

**7.3.** Using the pocket PC, create the list of accessions ready for thawing. Fill CF 30 Cool Box System with LN. Using sharp-pointed pliers, transfer viability assessment vials of 12 accessions from the small cryotank to CF 30 Cool Box System holder. Place cryovials in the holder in two columns of 6 accessions each. Separate columns with an empty column. Using sharp-pointed pliers, open the cryovials. Process accessions according their positions in the column.

**7.4.** Place 6 large filter papers (approx. 25x30 mm) on a stack of 3-6 sterile paper sheets. On each of the large filter papers, place 3 small sterile filter papers (approx. 10x10 mm). The "filter paper sandwich" absorbs the excess of RS after thawing. Place 3 sterile large filter papers (approx. 25x30 mm) on each petri dish with sterile Potato Recovery Medium I (PRM-0.3 M, see OP108). The large filter papers support the shoot tips during the first three stages of the recovery process (Recovery Medium I to III, see OP108). Label petri dishes with remaining RS label. The order of accessions within the stack of petri dishes should coincide with the order the accessions are processed (= position of accession within column of cool box).

**7.5.** Using cotton forceps, quickly transfer the 3 aluminum strips of the first accession from the cryovial to the corresponding RS test tubes (one strip per test tube). Remove empty cryovials from freezing container. Place the accession's cryovial label on the corresponding recovery cycle petri dish as an additional quality control (the label has to dry first).

**7.6.** Immediately shake the test tube to ensure rapid thawing. Take care that all shoot tips are submerged in RS and avoid sticking to the test tube wall. Treat shoot tips for 20 minutes in RS (at room temperature). Note the finishing time of thawing with RS on the sticker and monitor the RS exposure time with a digital clock.

## **8. Recovery Cycle**

**8.1.** Using the RS pipette, suck up the shoot tips from RS and place them for approx. 1-3 minutes on small sterile filter papers to drain excess of RS. For transfer of shoot tips from draining paper to recovery media, carefully unload shoot tips slipping the small filter papers, with shoot tips oriented downwards, over one of the large filter papers supported on Recovery Medium I of the respective accession. Check that labcode on petri dish label and the RS test tubes match.

**8.2.** Seal the petri dish with parafilm. Identify petri dish with original cryovial label. After completing 6 accessions, scan petri dish labels and record thawing data with the pocket PC. Using the hand-held barcode printer, print detailed labels for petri dish identification. Label petri dishes. Labels contain the following information: accession number, labcode, ID of person that cryopreserved the accession, date of freezing, date of thawing, and the date shoot tips were placed on Potato Recovery Medium I. Wrap petri dishes with aluminum foil (darkness). Incubate shoot tips for 3 days at  $20\pm 1$  °C in darkness.

**8.3.** After 3 days on Recovery Medium I, transfer filter papers holding shoot tips to sterile Recovery Medium II (PRM-0.2 M, see OP108). Transfer the lid of the original petri dish, the Recovery Medium I petri dish, to the new petri dish with Recovery Medium II (to facilitate labeling). Scan original petri dish labels with pocket PC, record transfer data, and print new labels with the hand-held barcode printer. Labels contain the following information: accession number, labcode, ID of person that cryopreserved the accession, date of freezing, date of thawing, and the dates shoot tips were placed on Potato Recovery Medium I and II. Place new label over original label. Wrap petri dish in aluminum foil (darkness) and incubate for 3 days at  $20\pm 1$  °C.

**8.4.** After 3 days on Recovery Medium II, transfer filter papers holding shoot tips to sterile Recovery Medium III (PRM-0.1 M, see OP108). Transfer the lid of the original petri dish, the Recovery Medium II petri dish, to the new petri dish with Recovery Medium III (to facilitate labeling). Scan original petri dish labels with the pocket PC, record transfer data, and print new labels with hand-held barcode printer. Labels contain the following information: accession number, labcode, ID of person that cryopreserved the

accession, date of freezing, date of thawing, and the dates shoot tips were placed on Potato Recovery Medium I, II and III. Place new label over original labels. Wrap petri dish with aluminum foil (darkness) and incubate for 3 days at  $20 \pm 1$  °C.

**8.5.** After 3 days on Recovery Medium III, transfer shoot tips to sterile Recovery Medium IV (PRM-25 g, see OP108). Before transferring an accession, scan the petri dish barcode label with the pocket PC and print 3 labels per accession (using the hand-held barcode printer). At this stage, place the shoot tips directly on the culture medium without a filter paper support. Transfer the 10 shoot tips supported on each of the filter papers to an individual petri dish with Recovery Medium IV. The labels contain the following information: accession number, labcode, ID of person that cryopreserved the accession, date of freezing, date of thawing, the dates shoot tips were placed on Potato Recovery Medium I, II and III and IV, IDs of the cryorun (1<sup>st</sup> or 2<sup>nd</sup>), and the petri dish number ("1", "2" or "3"). Incubate shoot tips for 4 days at  $20 \pm 1$  °C under diffuse light with a photoperiod of 16h/8h (light/darkness). For diffuse light, place a aluminium paper sheet onto the petri dishes. Petri dishes with signs of fungal/bacterial contamination are processed as described in step 10, without transferring the shoot tips.

**8.6.** After 4 days, remove the aluminum foil from the petri dishes and incubate shoot tips for approx. 18 additional days at  $20 \pm 2$ °C and  $85 \pm 20$   $\mu\text{mol.m}^{-2}\text{s}^{-1}$  light intensity, provided by cool daylight fluorescent tubes with a 16h/8h photoperiod.

**8.7.** 30-35 days after thawing, transfer the surviving and recovered shoot tips to 13x100 mm test tubes with Recovery Medium IV. Placing one sample per tube (see step 9).

## **9. Routine Viability Evaluations (30-35 and 60-65 days after thawing)**

**9.1.** Evaluate survival and re-growth rates 30-35 and 60-65 days after thawing. In the first evaluation, 30-35 days after thawing, a shoot tip is classified as survived when green tissue is observed, and as recovered when a shoot grows directly or laterally from the meristematic dome without any intermediate callus formation. Using a permanent marker, label survived shoot tips with a black unfilled circle on the base of the petri dish (externally). Label recovered shoot tips with a filled black circle. Dead or phenolized samples are not marked at all. Evaluate each petri dish separately. Scan the petri dish label with the pocket PC and record the following data: date of evaluation, ID of person that has done the evaluation, total number of shoot tips per petri dish, an the number of survived and regrown shoot tips per petri dish. The ID of the cryo-run (1st or 2nd), and the petri dish number ("1", "2" or "3") are included in the bi-dimensional barcode label of the Recovery Medium IV petri dishes and are therefore recorded automatically. For all 3 petri dishes, obtain complete data. Process petri dishes with signs of fungal /bacterial contamination as described in step 10, without transferring the shoot tips to test tubes.

**9.2.** For transfer to test tubes, process each petri dish individually, 1-5 days after the first evaluation. Stack petri dishes ordered by labcode of the accession. Using the pocket PC scan the barcode label of the petri dish that contains the shoot tips to be transferred. The display of the pocket PC shows the evaluation data, ID of the petri dish, and the run number. Check if physical identifiers (black circles) coincide with recorded evaluation data. Using the hand-held printer, print labels for 13x100 mm test tubes. Based on the evaluation data of 30-35 days , a database algorithm selects how many labels to print for each petri dish that is transferred. The labels contain the following information: accession

number, labcode, ID of person that performed the transfer, dates of freezing and thawing, current date, IDs of the cryorun (1st or 2nd) and the petri dish number ("1", "2" or "3"). Label test tubes before transferring samples.

**9.3.** Using long forceps, carefully transfer samples from petri dishes to 13x100 mm test tubes with Potato Recovery Medium IV (PRM-25 g, see OP108), placing 1 sample per test tube. Handle each petri dish separately and transfer only surviving and/or regrown shoot tips to test tubes. Rooted samples are transferred directly (= transplanting). If the sample has an elongated stem without rooting, slightly cut the sample's base before transferring to stimulate root formation. Directly transfer small but living samples (bud without stem elongation). Store test tubes in racks organized by labcode and petri dish number (in rows).

**9.4.** Perform the second and final viability assessment 60-65 days after thawing (in test tubes). Only samples that developed into complete, normal-looking *in vitro* plants with elongated stem, apex, leaves and roots are counted as recovered. For assessment, evaluate and record each of the different petri dish groups separately ("1", "2" or "3"). Using the pocket PC, scan the label of one test tube corresponding to a particular petri dish group. Check that the remaining test tubes correspond to the same accession, cryorun (1st or 2nd) and petri dish group ("1", "2" or "3"). Using the pocket PC record the following data: date of evaluation, ID of person that has done the evaluation, total number of shoot tips, and the number of survived and regrown shoot tips. The ID of the cryo-run (1st or 2nd) and petri dish group ID ("1", "2" or "3") are included in the bi-dimensional barcode label of the test tubes and are therefore recorded automatically. Additionally, record the following information: number of complete plants (recorded as recovered shoot tips), number of plants without roots (to be cut at its base and transferred to fresh culture medium), number of plants with deformation (recorded as survived shoot tips), phenolized samples (recorded as dead shoot tips), number of samples with small buds (to be transferred directly to fresh culture medium), number of white/dead samples (recorded as dead shoot tips). Groups of test tubes with signs of fungal/bacterial contamination are processed as described in step 10 ("Contamination Recheck").

**9.5.** Transfer samples that can't clearly be classified as recovered (e.g. survived shoot tips), or that still haven't developed into a complete *in vitro* plant (i.e. missing roots) to 13x100 mm test tubes containing fresh Recovery Medium IV (according to the criteria described in step 9.4.). Re-evaluate the plants 3 weeks later and update the evaluation data accordingly.

## **10. Contamination recheck**

**10.1.** Exceptionally, accessions might show signs of contamination in one of the following stages of the recovery cycle: (a) shoot tip transfer to petri dishes with Recovery Medium III or IV, (b) viability evaluations after 30-35 or 60-65 days. Using the pocket PC, scan the barcode label of the contaminated petri dish(es) or test tube(s) and record the following information: type of contamination (fungal or bacterial), date of record of contamination, and the ID of person that recorded the contamination. The ID of the cryo-run (1st or 2nd) and petri dish ("1", "2" or "3") are included in the bi-dimensional barcode label of the petri dish(es) or test tube(s) and are therefore recorded automatically. In order to obtain complete and reliable viability data for all three petri dishes, remove the required additional cryovials from the working or temporal storage tanks. Record cryovial removal from the tanks with the pocket PC to update stock and location of the accession. For contamination recheck, evaluate the viability of the removed cryovials as described in step 9 ("Routine Viability Evaluation").

**10.2.** If the recheck confirms the contamination, discard the respective run of the accession from the working or temporal storage tanks and cryopreserve again. If the recheck does not show any signs of contamination, record the evaluation data of the rechecked vials as original routine data, replacing the data of the contaminated vials. The removal of vials for contamination recheck can be done in parallel with the processes of discarding of accessions (see step 13) and transferring of accessions to the cryobank and temporary safety back-up tanks (see step 14).

### **11. Transferring of accessions to working tank (in cryo-laboratory)**

**11.1.** During the complete transfer process, personnel should wear ultra-low temperature gloves, cryo-apron, and safety goggles. For transferring, place a lab ice pan (capacity of 9 liter) in the laminar flow chamber and fill the container with approx. 3.5 cm of liquid nitrogen. If a transitory cryobox is labeled the first time (or re-labeled), collocate its cryobox ID on the four sides of the box (using ultra-low temperature resistant pencil). Using the pocket PC, print 4 cryo-resistant labels per cryobox and place them onto the hand-written identifiers (for security reasons, cryoboxes are identified both with cryogenic marker and bidimensional cryo-resistant barcode label). Unique IDs of transitory cryoboxes can range from Temp\_001 to Temp\_071 (for details, see next section on "Transfer of accessions to transitory tanks 1 and 2").

In general, the transitory cryo-boxes re-circulate between transitory tank 1 or 2 (Taylor Wharton LS-4800) and the cryo-box based working tank (CHART Cryosystem 2000). Consequently, they are already labeled adequately and can be used without additional identification.

Scan the cryobox label with the pocket PC. Place the cryobox in the lab ice pan with LN and submerge it half way into LN (using long forceps) without splashing, spilling or coming in contact with LN. If necessary, fill additional LN into the pan.

**11.2.** Place the cryo-cane based (CHART MVE XC47/11) and cryo-box-based working tank (CHART Cryosystem 2000) next to laminar flow chamber. The cryo-cane based tank contains the accessions to be transferred (see step 6.6.) and the cryo-box based tank will receive the accessions (stored in cryo-boxes). The cryobox-based tank has 4 box holders with a maximum capacity of 5 cryoboxes per holder (total capacity of tank: 20 cryoboxes). Place a maximum of 5 accessions in the transitory cryobox (12 cryovials per accession). The standard cryoboxes can hold 100 cryovials. Accessions occupy two rows within the cryobox. The remaining free positions within a row are not occupied by cryovials of another accession. For example: the first accession in the cryobox will occupy the positions 1 to 12, but the 8 free positions of the second row (position 13-20) will not be occupied. The second accession of the box will occupy the positions 21-32, and so on.

**11.3.** Using sharp-pointed pliers, quickly transfer accessions from the cryo-canes to the cryobox. Check that cryovials of the same accession have the same label. When the box is completely full, scan the first vial of each accession (positions 1, 21, 41, 61, and 81) with the pocket PC. The pocket PC will automatically assign the position to each accession (i.e. 1-12, 21-32, 41-52, 61-72 and 81-92). Check that the physical and electronically recorded numbers and positions of the cryovials coincide. Using the pocket PC, also record and check the following information: date of transfer, ID of the person who transferred the accession, ID of the cryorun of each accession (1st or 2nd), and number of cryopreserved shoot tips.

**11.4.** Carefully close the cryobox with its lid. A second person lifts a holder (with free positions) from the cryo-box based tank and removes the security bar, while the first person lifts the closed cryobox by one of its sides with long forceps, takes the box at its lifted side and places it in the corresponding slot of the tank holder. As soon as the box is in the slot, the second person inserts the security bar and re-submerges the holder in the tank. Close the tank. Proceed the same way with the other accessions to be transferred.

## **12. Transferring accessions from working tank (CHART Cryosystem 2000) to transitory tanks 1 or 2 (Taylor Wharton LS-4800)**

**12.1.** When the cryobox-based working tank is full (max. 20 cryoboxes), transfer boxes to the transitory tanks 1 and 2 (Taylor Wharton LS-4800). The transitory cryotanks have 6 holders each, with a capacity of 8 cryoboxes per holder. The total capacity per tank is 48 cryoboxes. The location of the cryoboxes within the transitory cryotanks is pre-defined. Transitory tank 1 holds the cryoboxes "Temp\_01" to "Temp\_48" and transitory tank 2 holds the cryoboxes "Temp\_49" to "Temp\_71". The transitory cryoboxes "Temp\_01" to "Temp\_08" are located in holder 1 of Transitory tank 1. The cryobox "Temp\_01" is located in the highest position of holder 1 (next to the tank's lid). "Temp\_02" in the second upper position, "Temp\_03" in the third upper position, and so on. The transitory cryoboxes "Temp\_09" to "Temp\_16" are located in holder 2 of transitory tank 1. "Temp\_09" is located in the highest position, "Temp\_10" in the second upper position, "Temp\_11" in the third upper position, and so on.

In this manner, the position of each of the transitory cryoboxes ("Temp\_01" to "Temp\_71") is exactly pre-defined, both within the transitory tanks and their holders. Holders 4 to 6 of transitory tank 2 are used for the transitory storage of sweetpotato accessions.

**12.2.** Place working tank (CHART Cryosystem 2000) next to transitory cryotanks 1 and 2 (Taylor Wharton LS-4800). While a second person lifts for a short moment the holder of the working tank, so that the first person can scan the barcode labels of the cryoboxes (using a pocket PC). An algorithm of the database automatically assigns the corresponding tank ID (T1 or T2) and holder (1 to 6) to the cryoboxes. Therefore, the boxes' position is pre-defined and fixed. Record and check the date of transfer and ID of person that performed the transfer with the pocket PC.

Transfer transitory cryoboxes from the working tank (CHART Cryosystem 2000) to their corresponding positions within the transitory cryotanks 1 and 2 (Taylor Wharton LS-4800).

## **13. Discarding accessions**

**13.1.** Accessions may be discarded from the cryotanks for different reasons: (a) confirmed fungal /bacterial contamination, (b) low routine recovery rates, (c) viability reassessment shows low recovery rate, (d) sub-optimal storing conditions, or (e) low stock.

**13.2.** Using the pocket PC, load the list of accessions to be discarded. Locate the tank and cryobox that contain the accession(s) to discard. Transfer the cryobox to a lab ice pan (capacity of 9 liters) that was previously filled with approx. 3.5 cm of LN. Scan the barcode label of the accession to be discarded (only the 1st cryovial of the corresponding accession). Using the pocket PC, record the date of discarding, the ID of the person who discarded the accession, and reason for discard. A record of all historical data will

be automatically saved in the database. Physically discard the corresponding cryovials. If no additional accessions will be discarded from the cryobox, close cryobox and place it back to its original location. Cryoboxes and tanks are handled as described before.

Discarding of accessions can be done in parallel with the processes of contamination recheck (see step 10) and transferring of accessions to the cryobank and temporary safety back-up tanks (see step 14).

#### **14. Transferring of accessions to cryobank and temporary safety back-up tanks**

**14.1.** Transfer accessions that meet the minimum recovery and quality criteria described in the step 1 from the transitory tanks (T1,T2) to the cryobank (P1, P2) and temporary safety back-up tank (BB-1). Boxes for the cryobank and safety back-up have unique identifiers that range from P001 to P480 and BBP001 to BBP180, respectively, where "P" stands for Potato and "BBP" for Black Box Potato (= Safety Back-up). Label cryobank and safety-back up boxes with cryogenic marker. Using the pocket PC, print 4 cryo-resistant labels per cryobox and place them onto the hand-written identifiers. For security reasons, cryoboxes are identified both with cryogenic marker and bidimensional cryo-resistant barcode label.

**14.2.** The cryobank and safety back-up tanks (CHART MVE-819P-190AF) have 4 sections. Each of the sections can hold 3 holders for standard cryoboxes (100 vials/box) and 1 holder for small cryoboxes (25 vials/box). Both types of holders have a capacity of 15 cryoboxes. The total capacity of one of the CHART MVE-819P-190AF tanks is 180 standard boxes (100 vials/box) and 60 small cryoboxes (25 vials/box). The location of the cryoboxes within the cryobank and safety back-up tanks are pre-defined, similar as described before for the Transitory Tanks 1 and 2, with the difference that the order of the cryo-boxes within the holder are inverted.

For example: The cryobank boxes "P001" to "P015" are located in holder 1 of cryobank tank 1 (P1), where "P001" is located in the lowest position (next to the bottom of the tank), "P002" in the second lowest position, "P003" in the third lowest position, and so on. The cryobank boxes "P016" to "P030" are located in holder 2 of cryobank tank 1 (P1), where "P016" is located in the lowest position (next to bottom of tank), "P017" in the second lowest position, "P018" in the third lowest position, and so on. The 100-vial cryobox holders have IDs ranging from 1-3, 5-7, 9-11 and 13-15, while the 25-vial box holders have the IDs 4, 8, 12 and 16. The temporal safety-back up tank (BB-1) is organized in a homologous manner, but cryoboxes are identified with "BB"-prefixes. In the safety-back up tank, the holders 1-3, 5-7 and 9-11 (100-vial boxes) and 4, 8, and 12 (25-vial-boxes) are assigned to potato, while the holders 13-16 are used for the sweetpotato cryo-backup. First fill the 100-vial boxes, then the 25-vial boxes.

**14.3.** Make sure that the cryobox-based working tank (CHART cryosystem 2000) is empty and can be used as temporary storing facility during the transfers. Set up two lab ice pans (capacity of 9 liter) and fill them with approx. 3.5 cm of LN. Place cryobank and safety-back up boxes (without lids) in the closest pan and submerge them half way into LN (using long forceps) without splashing, spilling or coming in contact with LN. If necessary, pour additional LN into the pans. Using the pocket PC, scan the cryobox ID labels of cryobank and safety-back up box (= destiny).

Using the pocket PC, load the list of accessions ready to be transferred and order them descending by their cryobox IDs (of the transitory cryoboxes). Place the transitory cryobox (for example: Temp\_03) that contains accessions to be transferred in the other lab ice pan. Using the pocket PC, scan the cryobox ID label of transitory cryobox (= source).

Using the pocket PC, scan cryovial label of the accession to be transferred. The pocket PC will automatically show the distribution of cryovials between cryobank and safety back-up cryoboxes and their position in the cryoboxes. Using long or cotton forceps, transfer accessions from the temporary storage box to the cryobank and safety back-up boxes. Confirm transfer in the pocket PC.

**14.4.** Once all cryovials have been transferred to the cryobank and safety-back up boxes, close the transitory cryobox and return it to its position in the transitory tank (T1 or T2) or use the box-based working tank (CHART Cryosystem 2000) for temporal storage and transfer transitory cryoboxes by groups to T1 and T2. Transfer next cryobox (with accessions to transfer) from the transitory tank to the polystyrene container with LN and proceed as described before. It is not necessary to scan the cryobank or safety-back box labels again (only if they are replaced by a new empty one). Every time the transitory cryobox is replaced by a new one, scan the cryobox ID barcode label of the transitory cryobox. Always monitor the LN level in the polystyrene boxes and refill if necessary.

**14.5.** When the cryobank or safety back-up box is full, check the following information before closing the box with its lid: date of transfer, person who performed the transfer, number of vials and samples per box, ID of cryobox, position of cryobox, position of vials within cryobox, and the cryorun ID (1st or 2nd). All this information is filled out automatically by an algorithm, but it is necessary to check that the electronic data, physical position, and identification coincide. Transfer the cryobank or safety back-up box to the box-based working tank (CHART Cryosystem 2000) for temporal storage, and transfer them at the end of the working day to the cryobank and safety back-up tanks (less manipulation of cryobank and safety back-up tanks). Continue processing all accessions that meet the criteria for transfer. Every time an accession is transferred, it is automatically highlighted in the transfer list of the pocket PC (additional control).

Transferring of accessions to the cryobank and temporary safety back-up tanks can be done in parallel with the processes of contamination recheck (see step 10) and discarding of accessions (see step 13).

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