Genebank: In vitro conservation of potato - OP025

Date : September 02nd, 2016

Reference : Panta, A.; Cruzado, J.; Silvestre, R.; Franco, N.; Barkley, N.A.; Ellis, D. 2015. Genebank: In vitro conservation of potato. CIP-OP025. v.181:15 pps. 1-26 pp.

Not for general distribution.

For most current version, please contact Ana Panta (a.panta@cgiar.org)

Genebank: In vitro conservation of potato, v.181. CIP-OP025 by International Potato Center (CIP) is licensed under a Creative Commons Attribution-ShareAlike 4.0 International License. Permissions beyond the scope of this license may be available at <http://www.cipotato.org/>

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| | |
|------------------------------------|---|
| TITLE | Genebank: In vitro conservation of potato - OP025 |
| OWNER * | In vitro Conservation Curator |
| APPROVER * | Head Genebank |
| APPROVAL DATE * | Sep 14, 2012 |
| LAST REVIEW DATE * | September 02nd, 2016 |
| REVIEW FREQUENCY * | 12 Months |
| ISSUE DATE | Sep 05, 2016 23:36 |
| CONTRIBUTORS * | Panta, Ana (CIP), Silvestre, Rocio (CIP), Franco, Nataly Elizabeth (CIP) |
| KEYWORDS | accredited, procedure, invitro-conservation |
| DOCUMENT ID * | OP025 |
| VERSION NUMBER | v.181 |
| COMPETENT PERSONNEL * | A.Panta; R.Silvestre; N.E.Franco; C.M.Ramirez; M.Yupa; J.Cruzado; J.Sanchez; M. Marin; A.Alvarado; V.Mallma;I.Espinoza; R.Serrato |
| RELATED CLAUSES | 5.8 |
| RELATED DOCUMENTS | RD017;RD013;RD029 |
| CITATION OF THIS DOCUMENT * | Panta, A.; Cruzado, J.; Silvestre, R.; Franco, N.; Barkley, N.A.; Ellis, D. 2015. Genebank: In vitro conservation of potato. CIP-OP025. v.181:15 pps. |
| CONTACT | For most current version, please contact Ana Panta (a.panta@cgiar.org) |

- All fields marked with * are required

ENGLISH VERSION

INTRODUCTION

In vitro potato germplasm conservation utilizes tissue culture techniques reducing plantlets' growth. This is done applying osmotic stress treatments and incubation at low temperature. The methods described below allow the medium-term conservation of about 8,700 potato genotypes, following subcultures approximately every two years.

SCOPE

The scope of this procedure covers the conservation of all *in vitro* potato germplasm material in the genebank.

SAFETY

Specific requirements above laboratory safety procedures are described in the document "Good practices during in vitro bank activities", item 2.3.

Each member of staff is obligated to read, learn and applied safety rules following the above document. It is mandatory that during labor inside the cold culture chamber ($7\pm 1^{\circ}\text{C}$) the staff wear warm protecting coat, and the labor must be realized in 30 min periods, after each 30 min period the staff must leaf from the chamber at least 15 min.

MATERIALS

| EQUIPMENT | |
|--------------------------------------|-----------------------|
| Autoclave | Medium dispenser |
| pHmeter | Analytical balance |
| Laminar flow chamber | Magnetic Stirrer |
| Refrigerator | Glass bead sterilizer |
| Pocket Personal Computer (Pocket PC) | Bar coding scanner |
| Desk Computer | Tools sterilizer |
| Freezer (-20°C) | Precision balance |

| EQUIPMENT | |
|--|---------------------------------|
| Water deionizer | Destillator |
| Microwave oven | Sterilizing oven |
| Environmental monitoring device (HOBO data-logger) | Incubator Chamber (20° ±2°C) |
| Thermal printer for bar-code labelling | Cold Chamber (7° ±1°C) |
| Drying oven | |
| OTHER MATERIALS | |
| Glass test tubes (25x150 mm, 16x125 mmm) | Cotton |
| Test tube caps | Labels |
| Forceps | Wash bottle with alcohol (70%) |
| Scalpels | Burner with alcohol (96%) |
| Blades No. 10 | Sterile paper sheets |
| Saran Wrap | Test tube racks |
| Alcohol Gel | Tool sterilizing box |
| Tool support | Glass vessel with alcohol (96%) |
| Culture media (see table 3) | Parafilm |

PROCEDURE

General rules:

- Conservation can be long- or transitory-term according to scientists' request and material type. Germplasm (G) and pathogen tested materials (PT), as well as some accessions of research material (R) , are conserved long-term. However most research material is conserved for a transitory-term. The name of the responsible scientist and the storage period for each accession have to be recorded in the Tissue Culture Laboratory database (TCL-DB).

- Accessions are grouped by health status: HS0, HS1, or HS2, based on the germplasm health status CIP's categorization (RD013); and HS2R, comprising materials HS2 requiring health status re-check. Each group comprises three material types:
 - (a) **'Germplasm'**, containing native cultivars, and wild genotypes from accessions requiring clonal propagation (non-producing botanical seed); these accessions are identified with the letter "G" in their corresponding laboratory code (Labcode).
 - (b) **'Pathogen tested'**, containing improved varieties and cultivars representing the crop diversity (core collection), made pathogen tested until 1999; they are identified by the letters "PT" in the Labcode.
 - (c) **'Research material'** comprises all material in trust not corresponding to any group described before; they are mainly breeding materials and are identified with the letter "R" in their Labcode.

1. Storage:

1.1. *In vitro* cultures are stored in a chamber under controlled conditions: temperature of $7 \pm 1^\circ\text{C}$, light intensity of $5\text{-}20 \mu\text{mol}\cdot\text{m}^{-2}\cdot\text{seg}^{-1}$ with a 16h/8h photoperiod of light and darkness (Fluorescent lamp COOL DAYLIGHT, 36W) (*). Under these conditions, *in vitro* cultures can be viable for approximately two years.

1.2. Three tubes per accession are conserved, with four plantlets per tube. Each tube must be labeled indicating the barcode and the following passport identification data: CIP's accession number, labcode, cultivar name, type and lot of culture medium, health status (HS0, HS1, or HS2), and subculture date.

1.3. Thirty accessions are stored in each test tube rack. Place three to four test tube racks per shelf section. Each test tube rack should be labeled indicating its location. Each shelf section has a number. The location of the accessions has to be recorded in the TCL-DB.

* Illuminance is constantly checked using an environmental monitoring device (HOBO logger). Fluorescents of low illuminance capacity, are replaced by new ones after approximately 3,000 h of use.

2. Evaluation of viability:

2.1. Evaluate the *in vitro* cultures every 3-4 months.

2.2. Transfer the *in vitro* cultures from the conservation chamber at $7 \pm 1^\circ\text{C}$ to the active collection chamber at $20 \pm 2^\circ\text{C}$. The test tube racks containing the accessions should be transported by laboratory cart.

2.3. Using a pocket PC record in the TCL-DB the accessions contained in each test tube rack (following the Data record and monitoring guide of the *in vitro* Genebank - point 3).

2.4. Verify that accessions location is recorded correctly in the TCL-DB. Verify that missing accessions and those with only one test tube in stock are in subculture process. The *in vitro* cultures of these accessions should be stored in the active collection chamber. Locate the accessions using a pocket PC (following the Data record and monitoring guide of the in vitro Genebank - point 2).

2.5. Evaluate the accessions being in subculture. This has to be done for:

- (a) Routine renewal
- (b) Vigor increasing treatments (see Annex 1)
- (c) Microbial contaminant elimination treatments (see Annex 1)

The evaluation is done visually according to the indicators specified in Table 1. When the *in vitro* cultures continue to have low viability and/or microbial contamination, repeat the specific treatments.

2.6. Continue evaluating all the other accessions of the test tube rack. Accessions classified as 'regular', as well as those of the 'medium' category with only two test tubes in storage, have to be sub-cultured for their renewal. For this purpose, one of the test tubes with 'medium' vigor status is transferred to a new test tube rack and placed in a transitory location. For all cases, using a pocket PC record the evaluation results, the new transitory location, and the next activity to follow (urgent renewal, renewal, conservation) in the TCL-DB (following the Data record and monitoring guide of the in vitro Genebank - point 2).

2.7. Return the test tube rack containing the evaluated accessions to the cold chamber. Using a pocket PC, verify that its location is correctly recorded in the TCL-DB.

2.8. Place the transitory test tube rack in the cold chamber in the corresponding site, according to the health status (HS0, HS1, HS2 or HS2R). Identify the test tube rack with a plate indicating the corresponding rack and shelf number.

2.9. The number of stocks and the location of all accessions contained in the transitory test tube rack are recorded in the TCL-DB.

2.10. Sub-culture the accessions in fresh medium no longer than one week after its evaluation.

3. Sub-culturing for renewal:

3.1. List and labeling:

3.1.1. Each test tube rack to process is assigned to an operator. Transfer the test tube rack, containing the accessions to be processed, to the transferring area.

3.1.2. Generate the list of accessions to be processed following the Data Monitoring Manual.

3.1.3. Record the culture medium to be used (i.e. MSA) in the TCL-DB.

3.1.4. Only generate labels for the accessions that can be processed during one day. DO NOT print additional labels for any reason. The operator has to process the accessions on the day the labels were printed. Print three labels per accession. Prepare labels following the Data record and monitoring guide of the in vitro Genebank - points 6 and 7

3.2. Sub-culture:

3.2.1. At least 30 minutes before starting the subculture, turn on and disinfect the laminar flow chamber, turn on the tools sterilizer, and prepare all required materials and tools following the protocol of Good practices during in vitro activities

3.2.2. Under aseptic conditions process accessions one by one; each accession in process is placed in one rack inside the laminar flow chamber. Use one tube per accession as original material for the subculture process. Each original tube contains four plantlets. Locate and separate the three labels corresponding to the specific accession to be sub-cultured.

3.2.3. Using a scalpel remove the tube's sealing tape (i.e. parafilm) and the label. Keep this old label beside the new labels of the accession. This label will be placed as a control on one of the new tubes.

3.2.4. Subculture the plantlets following the In vitro multiplication of potato - OP056 protocol. Use 16 x 125 mm test tubes with MSA medium (3 ml/tube) for the subculture. Four stem segments are placed into each test tube. Each stem segment contains 1-2 nodes. Subculture three tubes per accession.

3.2.5. Label the tubes. Verify that the old label (coming from the test tube that contained the mother plants) has the same identifications as the new labels.

3.2.6. Attach a saran wrap tape to each tube. Attach the labels on the saran wrap tape. The old label has to be attached to one of the three tubes, in the cap. Seal each tube around the cap with saran wrap tape.

3.2.7. After processing all accessions from one test tube rack, place the rack in its corresponding location in the active collection chamber. Incubation conditions are: temperature of $20 \pm 2^{\circ}\text{C}$, light intensity of $50\text{-}95 \mu\text{mol}\cdot\text{m}^{-2}\cdot\text{s}^{-1}$ with a 16h/8h photoperiod of light and darkness (Fluorescent lamp COOL DAYLIGHT, 36W) (*)

3.2.8. Asepsia monitoring:

(a) 5-7 days after the sub-culture, evaluate the *in vitro* cultures visually for verifying the absence of microbial contamination signs (fungi, bacteria).

(b) Using the Computer generate the accession list. Evaluate visually each accession and record the evaluation results in the TCL-DB (using the pocket PC), following the Data record and monitoring guide of the in vitro Genebank- points 5.1 and 5.2

(c) Contaminated cultures are removed from the test tube rack and placed inside a bag for autoclaving. The bag is transferred to the washing area for its sterilization.

(d) If all the tubes from one accession are contaminated, repeat the subculture process using a new stock tube from the conservation chamber. When all stocks are showing bacteria contamination, the accession must be transferred to the phytosanitary area where will be processed for bacteria elimination.

(e) Accessions with no signs of contamination will continue their incubation in the active collection chamber.

3.2.9. Viability monitoring:

(a) 20-25 after the subculture evaluate the viability of the *in vitro* cultures according to the indicators of Table 2. Using the Computer generate the accession list. Evaluate visually each accession and record the evaluation results in the TCL-DB (using the pocket PC), following the Data record and monitoring guide of the in vitro Genebank - points 5.1. and 5.3

(b) Accessions classified as 'medium' or 'good' are transferred to a new test tube rack and placed in a transitory location in the chamber $20^{\circ}\pm 2^{\circ}\text{C}$. These accessions are ready for next subculture in MSA (when more material is required for safety copies) or on conservation medium (MSA + Sorbitol).

3.2.10. Using the pocket PC record in the TCL-DB the location of the test tube rack which contains the accessions that are ready for its sub-culture on conservation medium; and should be processed as soon as possible.

* Luminance is recorded by a environment conditions data logger (HOBO). Fluorescents of low luminance capacity are replaced by new ones after approximately 3,000 h of use.

3.3. Transfer to conservation medium:

3.3.1. Accessions that are ready for its subculture on conservation medium were placed in test tube racks in the active collection chamber (see point 3.2.1 to 3.2.10). Each test tube rack to be processed is assigned to an operator. Move the test tube rack, containing the accessions to be processed, from the active collection chamber to the transferring area. Each rack will contain only 5 accession; accessions are placed on rack cells keeping 2 empty rows between each accession, this practice minimize risk or errors in mixing tubes.

3.3.2. Generate the accession list following the manual of data monitoring.

3.3.3. Record the medium change, ie. from MSA to the corresponding conservation media (S42, S32 o S22), in the TCL-DB following the Data record and monitoring guide of the in vitro Genebank - point 4

3.3.4. The conservation medium to be utilized is defined according to the viability of the accession on MSA medium. The utilized culture media contain different concentrations of sorbitol, which is a growth retardant (see Table 3). Accessions can be subcultured in three ways:

(a) Accessions evaluated as 'good': all tubes are transferred to S-42 medium.

(b) Accessions evaluated as 'medium': all tubes are transferred to S-32.

(c) Accessions of wild species evaluated as 'medium': all tubes are transferred to S-22.

3.3.5. Only generate labels for the test tube rack that is processed. DO NOT print additional labels for any reason. The operator has to process the accessions on the day the labels were printed. Print three labels per accession. Prepare labels following Data record and monitoring guide of the in vitro Genebank - points 6 and 7

3.3.6. Under aseptic conditions process accessions one by one; each accession in process is placed in one rack inside the laminar flow chamber. Separate the labels corresponding to the accession to be sub-cultured in conservation medium. Using a scalpel remove the tube's sealing tape (i.e. saran wrap) and the label. Keep the old labels (subculture on MSA medium) beside the new labels of the accession. These labels will be placed as control on the new tubes, in the caps.

3.3.7. Subculture the plantlets following the *In vitro* multiplication of potato - OP056 protocol. Use 25 x 150 mm test tubes with S-22, S-32 or S-42 medium (10ml / tube) according to the evaluation results (see point 3.3.4.). Four apical stem segments are placed into each test tube. Each apical stem segment contains 1-2 nodes. Subculture three tubes per accession. Always prefer apical buds; exceptionally, axillary buds can be used.

3.3.8. Label the tubes. In order to avoid any mislabeling risk, it is mandatory to verify that the all old labels (coming from the tubes containing the mother cultures) has the same identifications as the new labels.

3.3.9. Attach a Parafilm tape to each tube. Attach the labels on the saran wrap tape. The label from one of the tubes that contained the mother cultures must be attached to one of the new tubes; this label is used as identification witness. Seal the tubes with Saran wrap tape.

3.3.10. Asepsia monitoring:

(a) 5-7 days after the subculture, evaluate the *in vitro* cultures visually for verifying the absence of microbial contamination signs (fungi, bacteria).

(b) Using the Computer generate the accession list. Evaluate visually each accession and record the evaluation results in the Genebank-database (using the pocket PC), following the Data record and monitoring guide of the in vitro Genebank - points 5.1 and 5.2

(c) Contaminated *in vitro* cultures are removed from the test tube rack and placed inside a bag for autoclaving. The bag is transferred to the washing area for its sterilization.

(d) If all the tubes from one accession are contaminated, repeat the subculture process using a new stock tube from the conservation chamber.

(e) Accessions with no signs of contamination will continue their incubation in the active collection chamber.

3.3.11. Viability monitoring:

(a) 15-20 days after the subculture evaluate the viability indicators according to Table 2. Using the Computer generate the accession list. Evaluate visually each accession and record the evaluation results in the TCL-DB (using the pocket PC), following the Data record and monitoring guide of the in vitro Genebank- points 5.1 and 5.3

(b) Accessions classified as 'medium' or 'good' are transferred to a new test tube rack and placed in a transitory location in the active collection chamber. These accessions are ready for its transfer to the cold chamber.

3.3.12. Using the pocket PC, record in the TCL-DB the accessions ready to be transfered from the active collection chamber to the conservation chamber.

(a) Using the pocket PC, select the option "Save" (down site)

(b) Select the option "Move accessions to CF" (cold chamber)

3.3.13. Transfer the accessions to the conservation chamber and place them in its corresponding location. Using the PC, verify that the transfer date, stock numbers, and locations have been correctly recorded in the TCL-DB.

3.3.14. Evaluate the viability every 3-4 months following points 2.1 to 2.10) see Data record and monitoring guide of the in vitro Genebank - points 5.1 and 5.4. Additionally, evaluate the cultures monthly for verifying absence of vessels with microbial contamination.

Table 1. Indicators utilized in the visual viability evaluation of potato accessions for medium-term conservation (culture media S22, S32 or S42)

| Viability category | Stem necrosis | Chlorosis | Hyper-hydration | Etiolation |
|--------------------|---------------|-------------|-----------------|-----------------|
| Good | 0-10% | No | No | No |
| Medium | 10-30% | Slight | No | Slight |
| Regular 1 | 30-70% | No, slight | No-slight | No |
| Regular 2 | 0-30% | Slight-high | Slight-high | No, Slight-high |
| Dead | 100% | - | - | - |

Table 2. Categories utilized in the visual viability evaluation of sub-cultured potato accessions (culture medium MSA)

| Viability category* | *Chlorosis * | *Rooting * | Necrosis | Hyper-hydration | Shoot growing |
|---------------------|--------------|------------|----------|-----------------|---------------|
| Good | No | Good | 0-10% | No | Good |
| Medium | Slight | Medium | 10-30% | No-slight | Medium |
| Regular 1 | Slight-high | Poor | | Slight-high | Poor |
| Dead | - | - | 100% | - | - |

Table 3. Conservation media composition for potato *in vitro* culture (S22, S32, S42)

| Components | S22 | S32 | S42 |
|-----------------------|------|------|------|
| MS salts* (g/l) | 4.33 | 4.33 | 4.33 |
| Glycine-HCl (mg/l) | 2 | 2 | 2 |
| myo-inositol (mg/l) | 100 | 100 | 100 |
| Nicotinic acid (mg/l) | 0.5 | 0.5 | 0.5 |
| Pyridoxine-HCl (mg/l) | 0.5 | 0.5 | 0.5 |
| Thiamine-HCl (mg/l) | 0.1 | 0.1 | 0.1 |
| Sucrose (g/l) | 20 | 20 | 20 |
| Sorbitol (g/l) | 20 | 30 | 40 |
| Agar (g/l) | 6.5 | 6.5 | 6.5 |

| | | | |
|----|-----|-----|-----|
| pH | 5.6 | 5.6 | 5.6 |
|----|-----|-----|-----|

Murashige and Skoog basal Salts (1962). Provider: CAISSON.

INTERNAL QUALITY CONTROL

1. Test tubes are identified with labels with bar and bi-dimensional coding .
2. Racks containing *in vitro* cultures are identified with labels with bar coding.
3. *In vitro* material is evaluated for verifying absence of bacteria and fungi contamination.
4. Environmental conditions, temperature and light intensity, are monitored using Data logger (HOBO), and temperature is monitored using a remote system (SITRAD) .
5. For the stocks renewal of each accession, one label from the former *in vitro* cultures is attached to one of the new cultures.
6. There is a CIP's internal alarm system for supporting security of the personnel who works inside cold ($7\pm 1^{\circ}\text{C}$) culture chamber.
7. Internal audits and workshops on techniques and procedures revising are performed at least once per year.
8. Procedures are applied following validated methods and international standards (<http://www.biodiversityinternational.org/publications/search/results.html>).
9. Before starting the culture media preparation, verify that all reagents have not expired.

REFERENCES

1. Benson EE, Harding K, Debouck D, Dumet D, Escobar R, Mafla G, Panis B, Panta A, Tay D, Van den houwe I, Roux N. 2011. Refinement and standardization of storage procedures for clonal crops - Global Public Goods Phase 2: Part I. Project landscape and general status of clonal crop in vitro conservation technologies. System-wide Genetic Resources Programme, Rome, Italy.

<http://www.biodiversityinternational.org/publications/search/results.html>

2. Benson EE, Harding K, Debouck D, Dumet D, Escobar R, Mafla G, Panis B, Panta A, Tay D, Van denhouwe I, Roux N. 2011. Refinement and standardization of storage procedures for clonal crops Global Public Goods Phase 2: Part II. Status of in vitro conservation technologies for: Andean root and tuber crops, cassava, Musa, potato, sweetpotato and yam. System-wide Genetic Resources Programme, Rome, Italy.

<http://www.biodiversityinternational.org/publications/search/results.html>

3. Benson EE, Harding K, Debouck D, Dumet D, Escobar R, Mafla G, Panis B, Panta A, Tay D, Van den houwe I, Roux N. 2011. Refinement and standardization of storage procedures for clonal crops - Global Public Goods Phase 2: Part III. Multi-crop guidelines for developing in vitro conservation best practices for clonal crops. System-wide Genetic Resources Programme, Rome, Italy.

<http://www.biodiversityinternational.org/publications/search/results.html>

ANNEX 1

Treatments applied to increase vigour and to eliminate microbial contaminants

For all the accessions that require specific treatments transfer the stock with best healthy status to a transitory tube rack. The treatments are:

(a) Accessions showing poor growth and rooting are transferred to highly nutritive media (i.e. MPR, see table 4). The treatments applied to increase the vigour should be done as soon as possible.

(b) Accessions showing fungi with or without bacterial contaminants are disinfected using calcium hypochloride (2.5%) and transferred to fresh MSA medium.

(c) Accessions showing bacterial contaminants are transferred to MSA medium containing the antibiotic Cefotaxime (200 mg/l).

(d) Accessions showing bacterial contaminants resistant to Cefotaxime are disinfected with sodium hypochloride (2.5%) and transferred to MSA medium containing the antibiotic Rocephin (200 mg/l).

Processing of accessions with all stocks showing bacteria contamination are transferred to specific area dedicated to bacteria diagnosis, antibiotic treatment and cleaning. The processing steps "c" and "d" are performed in the phytosanitary area. The treatments for the elimination of microbial contaminants (fungal) should be applied in 1-2 days after detection.

Annex Table 1. Components of potato rapid propagation medium (MPR)

| MPR Components | Quantity |
|-----------------------------|-----------------|
| MS salts (g/l) * | 4.33 |
| Gibberellic acid (mg/l) | 0.1 |
| Glycine-HCl (mg/l) | 2 |
| Myo-inositol (mg/l) | 100 |
| Nicotinic acid (mg/l) | 0.5 |
| Pyridoxine-HCl (mg/l) | 0.5 |
| Thiamine-HCl (mg/l) | 0.1 |
| L-arginine (mg/l) | 4 |
| Calcium panthotenate (mg/l) | 2 |
| Putrescine HCl (mg/l) | 10 |
| Sucrose (g/l) | 30 |
| Coconut milk (ml/l) | 10 |
| Phytigel (g/l) | 3 |
| pH | 5.6 |

- Murashige and Skoog salts, 1962. Provider: CAISSON

ANNEX 2

Glossary of keywords:

*- *Culture medium: Substrate containing mainly the macro and micro-nutrients, vitamins and carbon sources (sugars) that are required for the optimal growth and development of any explant.

*- *Explant: Segment of piece from a plant organism. For example, protoplast, cell, tissue, organ, stem or leave segment.

- *In vitro* tissue culture: Techniques whereby an explant is grown aseptically in a chemically defined medium and incubated under controlled environmental conditions.

- Sub-culture: Process of *in vitro* tissue culture applied for the renewal of the *in vitro* cultures; explants from existing cultures (mother) are cultivated in fresh culture medium.

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REFERENCIAS

1. Benson EE, Harding K, Debouck D, Dumet D, Escobar R, Mafla G, Panis B, Panta A, Tay D, Van den houwe I, Roux N. 2011. Refinement and standardization of storage procedures for clonal crops - Global Public Goods Phase 2: Part I. Project landscape and general status of clonal crop *in vitro* conservation technologies. System-wide Genetic Resources Programme, Rome, Italy. <http://www.biodiversityinternational.org/publications/search/results.html>
2. Benson EE, Harding K, Debouck D, Dumet D, Escobar R, Mafla G, Panis B, Panta A, Tay D, Van denhouwe I, Roux N. 2011. Refinement and standardization of storage procedures for clonal crops - Global Public Goods Phase 2: Part II. Status of *in vitro* conservation technologies for: Andean root and tuber crops, cassava, Musa, potato, sweetpotato and yam. System-wide Genetic Resources Programme, Rome, Italy. <http://www.biodiversityinternational.org/publications/search/results.html>
3. Benson EE, Harding K, Debouck D, Dumet D, Escobar R, Mafla G, Panis B, Panta A, Tay D, Van den houwe I, Roux N. 2011. Refinement and standardization of storage procedures for clonal crops - Global Public Goods Phase 2: Part III. Multi-crop guidelines for developing *in vitro* conservation best practices for clonal crops. System-wide Genetic Resources Programme, Rome, Italy. <http://www.biodiversityinternational.org/publications/search/results.html>

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